Edwards

ThruPort Systems

INTRACLUDER INTRA-AORTIC OCCLUSION DEVICE

Instructions for Use

Directory

English (EN)	1
Français (FR) French	3
Deutsch (DE) German	5
Español (ES) Spanish	8
Italiano (IT) Italian	10
Nederlands (NL) Dutch	12
Dansk (DA) Danish	15
Svenska (SV) Swedish	17
Ελληνικά (EL) Greek	19
Português (PT) Portuguese	21
Český (CS) Czech	24
Magyar (HU) Hungarian	26
Polski (PL) Polish	28
Slovensky (SK) Slovak	30
Norsk (NO) Norwegian	33
Suomi (FI) Finnish	35
Български (BG) Bulgarian	37
Română (RO) Romanian	39
Eesti (ET) Estonian	42
Lietuvių (LT) Lithuanian	44
Latviešu (LV) Latvian	46
Türkçe (TR) Turkish	48

English (EN)

IntraClude Intra-Aortic Occlusion Device

Instructions for Use

CAUTION: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician.

See Figures at the end of the Instructions for Use.

DESCRIPTION

The IntraClude intra-aortic occlusion device (model number ICF100) is a 10.5 Fr (3.5 mm), triple-lumen, 100 cm long catheter with an elastomeric balloon near its distal tip designed to occlude the ascending aorta in order to partition the aortic root from arterial circulation. The balloon expands to occlude a range of aorta sizes from 20 to 40 mm. This device is designed to be used in the femoral approach with the Edwards EndoReturn (ER21B or ER23B) arterial cannula or the Edwards Introducer Sheath (IS19A). The shaft is provided with an extended strain relief, which tapers from 10.5 Fr to the remaining 9 Fr catheter, and is designed to prevent kinking at the hub and to avoid compressing the shaft when the hemostasis valve of the Edwards EndoReturn (ER21B or ER23B) arterial cannula or Edwards Introducer Sheath (IS19A) is firmly closed.

The large central lumen of the IntraClude device serves three functions: accommodating the guidewire, delivering cardioplegia solution to the aortic root, and venting fluid and air from the aortic root. The two remaining lumens serve as conduits for balloon inflation/deflation and aortic root pressure monitoring. The hub has two flexible extension tubes with an integrated luer connection to provide access for the accessories. The shaft is provided with markers to indicate the insertion depth. A blue Clamp-Lock device, provided on the extended strain relief portion, allows the IntraClude device to be locked in position.

The IntraClude device is packaged with the following accessories that function with the IntraClude device or provide connection to other equipment:

- A rotating hemostasis valve (RHV), and attached stopcock with green handle, connected to the large central lumen for guidewire insertion, fluid injection (i.e. cardioplegia delivery), and aortic root venting through the IntraClude device.
- A Y-connector, connected to the green stopcock on the RVH, with three tubing clamps, which enables alternation between cardioplegia solution delivery and aortic root venting.
- A 3-way stopcock with blue handle connected to the inflation/deflation extension of the IntraClude device for controlling the fluid pathway to the balloon inflation/deflation lumen.
- An 8 ft. (244 cm) pressure line, with a male and female luer and a blue line/marking, for balloon pressure monitoring.
- A 3-way stopcock with red handle connected to the aortic root pressure extension of the IntraClude device for controlling the fluid pathway to the aortic root pressure lumen.
- An 8 ft. (244 cm) pressure line, with a male and female luer and a red line/marking, for aortic root pressure monitoring.
- A 35 ml syringe for inflating the balloon and priming the IntraClude device.
- Two white luer caps for capping the balloon and aortic root pressure lines.
- A 200 cm, 0.038" (0.97 mm) J-tip guidewire to guide the IntraClude device through the vasculature.

HOW SUPPLIED

STERILE and NONPYROGENIC in an unopened and undamaged package. Do not use if the IntraClude device shows signs of damage (i.e., cuts, kinks, crushed areas), or if package is damaged or open.

INDICATIONS FOR USE

The IntraClude intra-aortic occlusion device is indicated for use in patients undergoing cardiopulmonary bypass. The IntraClude intra-aortic occlusion device occludes and vents the ascending aorta when the balloon is inflated. The device's central lumen allows delivery of cardioplegia to arrest the heart. The pressure lumen allows monitoring of the aortic root pressure.

CONTRAINdications

RELATIVE

- Moderate to severe peripheral or aortic atherosclerosis
- History of thoracic trauma

ABSOLUTE

- Aneurysm of the ascending aorta
- Severe aortic regurgitation

COMPLICATIONS

The following complications may occur during or following use of the IntraClude device:

- Injury to the aorta or artery through which it is inserted, including perforation and dissection
- Plaque embolization
- Damage to the aortic valve
- Stroke
- Allergic reaction to contrast medium
- Wound infection/sepsis
- Pain at the insertion site
- Hematoma at the access site
- Arrhythmia
- Arterial thrombosis
- Myocardial injury
- Cardiac failure
- Peripheral nerve damage
- Lower limb ischemia

WARNINGS AND PRECAUTIONS

This device is designed, intended, and distributed for single use only. Do not re-sterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Ensure proper levels of anticoagulant therapy are maintained prior to insertion of the cannula and throughout cardiopulmonary bypass, to reduce the risk of complications due to thrombus formation on or within the cannula, and in the blood stream.

Do not fill the IntraClude device balloon with air as embolism may occur. Always use the approved guidewire (packaged with the IntraClude device or sold separately). Failure to do so may cause device failure.

When placed in the body, the device should be manipulated only when the balloon is deflated and while under Transesophageal Echocardiographic (TEE) and/or fluoroscopic observation. For repositioning, the balloon should be partially deflated. Failure to do so may cause damage to the vessel.

If increased resistance is felt at any time upon insertion or removal of the guidewire, dilators, introducer, or cannula, stop and investigate the cause before continuing. Inability to easily advance or remove these devices may indicate vascular disease or injury. Closely examine the device position within the vessel using fluoroscopy and/or transesophageal echocardiography (TEE) prior to proceeding. Failure to do so may increase the risk of damage to the vessel and potential dissection.

Products are known to contain phthalates, which may be found in device materials containing plasticizers such as DEHP and BBP. High exposure to such phthalates during medical treatments in children and pregnant or nursing women may raise a concern. A review of available data and literature supports the conclusion that the benefits outweigh the overall residual risk.

Care should be used when handling the devices. Damage may result from kinking or stretching the devices. Do not use the device if kinked as this may result in inability to deflate the balloon or other device failures.

Proper surgical procedures and techniques are the responsibility of the medical profession. Described procedures are provided for informational purposes only. Each physician must determine the appropriate use of this device for each patient based on medical training, experience, the type of procedure employed, and the benefits and risks associated with device use.

INSTRUCTIONS FOR USE

REQUIRED EQUIPMENT

The IntraClude device should only be used in an operating room equipped for conventional cardiac surgery. Additional equipment should include:

- Transesophageal Echocardiography (TEE) and/or Fluoroscopy capabilities
- Cardiopulmonary bypass pump capable of delivering cardioplegia solution and venting the aortic root
- If using Fluoroscopy, radiopaque contrast medium, diluted with sterile physiologic solution at a 6:1 saline:contrast ratio, can be used.
- Aortic root pressure monitoring system capable of measuring pressure of at least 200 mmHg

- Balloon pressure monitoring system capable of measuring pressure of at least 500 mmHg

- Bilateral radial/brachial artery monitoring

WARNING: Failure to dilute the contrast media as recommended may result in formation of crystals of contrast medium that could delay or prevent balloon deflation.

HANDLING AND PREPARATION

Prior to opening the package, prepare the IntraClude device placement monitoring systems (TEE and/or fluoroscope). Examine the following:

- Descending thoracic aorta for plaque
- Ascending aorta for plaque
- The diameter of the aorta just above the Sinotubular Junction
- Aortic valve function

WARNING: The presence of aortic insufficiency may compromise the delivery of cardioplegia solution and may result in prolonged time to arrest the heart.

Follow sterile techniques during IntraClude device preparation and use.

PREPARE THE INTRACLUDER DEVICE FOR USE

Visually examine the IntraClude device and sterile packaging for evidence of damage (e.g. rips, tears, etc.). If damaged, DO NOT USE.

When shipped, all lumens, including the balloon lumen, contain air. This air must be displaced prior to insertion to ensure that only fluids fill the lumens. To displace air, follow the steps below:

1. Remove all components from the package and place on sterile prep table.
2. Visually inspect all the components for damage. Do not use if there are signs of damage.
3. Make sure that all pre-assembled luer connections are tight.
4. Prime the balloon and lumen:
 - a) Fill the syringe with 10 mL of sterile physiologic solution (making sure there is no air in the syringe) and connect to the open port of the blue stopcock. Hold the syringe vertically (with the syringe tip pointing down), and in a slow and controlled manner, pull back on the syringe plunger to remove air from the balloon and lumen.
 - b) While keeping the syringe connected and held vertically, in a slow and controlled manner, depress the syringe plunger to inject sterile physiologic solution into the balloon.
 - c) While keeping the syringe connected and held vertically, pull back on the syringe to remove air from the balloon. Note - the tip of the IntraClude device should be pointing downwards to facilitate air removal through the de-airing hole, which is located at the proximal end of the balloon.
 - d) While keeping the syringe connected and held vertically, depress the syringe plunger to inflate the balloon (making sure not to inject air into the balloon). Inspect for any air in the balloon. If there is air in the balloon, repeat the previous two steps until all air has been removed. Once all air has been removed, continue to the next step.
 - e) Pull back on the syringe and turn off the stopcock to the IntraClude device. The balloon and lumen are now primed.

WARNING: Ensure all air is removed. Point balloon tip downwards, to facilitate air removal. Air left in any lumen may embolize during or after insertion.

WARNING: The injection and depression of the syringe plunger should be performed in a slow and controlled manner so as to not entrain air in the IntraClude device. Air left in any lumen may embolize during or after insertion.

5. Prime the Root Pressure lumen:
 - a) Fill the syringe with sterile physiologic solution and connect to the open port of the red stopcock.
 - b) Flush the lumen with sterile physiologic solution until a continuous stream of sterile physiologic solution exits the tip. Turn off the stopcock to the IntraClude device.

6. Insert the J-tip of the guidewire into the Rotating Hemostasis valve (RHV), using the straightener provided with the guidewire to straighten the J-tip for ease of insertion into the RHV, and feed the guidewire into the IntraClude device until the tip of the guidewire is flush with the tip of the IntraClude device.

CAUTION: Use only the 0.038 in. (0.97 mm), 200 cm, J-tip guidewire with the IntraClude device (packaged with the IntraClude device or sold separately). Failure to do so could result in damage to the device and difficulty in placing the device in its proper position in the ascending aorta or patient injury.

7. Prime the cardioplegia lumen and Y-connector:
 - a) Remove the white cap from the side port on the green stopcock and turn off the stopcock to the Y-Connector.
 - b) Connect the syringe, filled with sterile physiologic solution, to the open side port on the green stopcock.
 - c) While holding the RHV down, flush and tap the RHV until no air bubbles are visible and a continuous stream of sterile physiologic solution exits the RHV. Tighten the RHV while flushing.
 - d) Continue flushing until a continuous stream of sterile physiologic solution exits the tip of the IntraClude device. Turn off the green stopcock.
 - e) Making sure that all clamps on the Y-connector are open, flush the Y-connector while holding the ends vertical. Close all tubing clamps once all air has been removed.
 - f) Remove the syringe and replace the white cap on the stopcock.

WARNING: Ensure all air is removed. Point balloon tip downwards, to facilitate air removal. Air left in any lumen may embolize during or after insertion.

WARNING: Do not over tighten the RHV in order to prevent damage to the guidewire during insertion.

8. Prime the balloon pressure monitoring line:
 - a) Connect the 244 cm (8 ft.) pressure line with the blue marking to the blue stopcock.

Released Date: 2018-04-23

Printed Date: 2018-12-05

Status = Released

- b) Remove the white cap from the side port on the blue stopcock and connect the syringe, filled with sterile physiologic solution, and flush the pressure line until a continuous stream of sterile physiologic solution exits the end.
- c) Cap the end of the pressure line using the supplied cap.
- d) Replace the cap on the stopcock.

9. Prime the aortic root pressure monitoring line:

- a) Connect the 244 cm (8 ft.) pressure line with the red marking to the red stopcock.
- b) Remove the white cap from the side port on the red stopcock and connect the syringe, filled with sterile physiologic solution, and flush the pressure line until a continuous stream of sterile physiologic solution exits the end.
- c) Cap the end of the pressure line using the supplied cap.
- d) Replace the cap on the stopcock.

WARNING: Prior to initiating cardiopulmonary bypass, completely expel air from system. Ensure that all stopcocks are closed to prevent air from entering the system. Failure to eliminate air from the system could result in an air embolus.

INTRACLUDER DEVICE INSERTION AND DELIVERY

Once the IntraClude device has been properly prepared, it can be inserted into the patient. The IntraClude device is to be used with the Edwards EndoReturn arterial cannula (ER23B or ER21B) or the Edwards Introducer Sheath (IS19A).

1. Heparinize the patient and place the Edwards EndoReturn arterial cannula or the Edwards Introducer Sheath according to its Instructions for Use.
2. Using a tubing clamp, clamp the non-reinforced section of tubing attached to the hemostasis valve of the cannula. This isolates the hemostasis valve from arterial pressure.
3. Wet the IntraClude device balloon surface with sterile physiologic solution to facilitate introduction through the hemostasis valve.

WARNING: Check the balloon prior to insertion to ensure that a vacuum has been maintained. A loss of balloon vacuum indicates that there is a leak in the system which may result in unexpected balloon deflation and loss of aortic occlusion. If a leak is detected, replace the IntraClude device or leaking component prior to use.

4. Fully open the arterial cannula hemostasis valve and insert the tip of the IntraClude device into the cannula until the tip is against the tubing clamp. Ensure that the balloon has cleared the seal of the cannula hemostasis valve.

CAUTION: Always wet the balloon surface with sterile physiologic solution prior to insertion to facilitate IntraClude device introduction through the hemostasis valve or damage to the IntraClude device may occur.

CAUTION: Ensure that the cannula hemostasis valve is fully opened before inserting the IntraClude device or damage to the IntraClude device may occur.

5. De-air the arterial cannula by momentarily releasing the tubing clamp, allowing for back bleeding.

WARNING: Prior to initiating cardiopulmonary bypass, completely expel air from the system. Failure to eliminate air from the system could result in an air embolus.

6. Tighten the arterial cannula hemostasis valve to allow easy advancement of the IntraClude device but without leakage of blood.

CAUTION: Do not over tighten the cannula hemostasis valve in order to prevent crushing of the IntraClude device. If there is difficulty in advancing the IntraClude device, loosen the hemostasis valve.

7. Remove the tubing clamp from the cannula. Ensure the low flow marker (Figure 1) is visible outside of the hemostasis valve. This indicates that the balloon has not yet entered the perfusion flow path of the cannula.

8. While keeping the IntraClude device in place, advance the guidewire up to the descending aorta and verify that the guidewire is intraluminal through TEE and/or Fluoroscopy visualization. Continue to advance the guidewire up to the ascending aorta using TEE and/or Fluoroscopy visualization.

WARNING: Do not advance guidewire if resistance is felt. Inability to easily advance guidewire may indicate vascular disease or injury.

WARNING: If the guidewire is not visible in the ascending aorta it may have entered into the arch vessels. Retract the guidewire into the descending aorta and re-advance it to the ascending aorta.

WARNING: Do not advance the guidewire past the aortic valve as damage to the valve may occur.

9. While keeping the guidewire in place, bring the IntraClude device up to the descending aorta and verify that it is intraluminal through TEE and/or Fluoroscopy visualization. While keeping the guidewire in place, bring the IntraClude device up to the ascending aorta using TEE and/or Fluoroscopy visualization. The balloon should be located just above the sinotubular junction.

WARNING: Inadvertent advancement of the guidewire can occur during placement of the IntraClude device. Correct position of the guidewire should be monitored and repositioned as required to avoid injury to the aorta and surrounding structures.

WARNING: Ensure the guidewire precedes the IntraClude device to prevent damage to aortic valve or artery. Failure to advance over a guidewire may result in dissection or perforation of the artery. If a dissection or perforation is detected, do not initiate or continue bypass as it may propagate dissection.

WARNING: Do not advance the IntraClude device if resistance is felt. Inability to easily advance the IntraClude device may indicate vascular disease or injury.

WARNING: Do not advance the IntraClude device past the aortic valve as damage to the valve may occur.

10. Fully remove the guidewire from the IntraClude device. Allow back bleed from the IntraClude device RHV to remove any air. Once all the air has been removed, close the RHV on the IntraClude device.

WARNING: Failure to fully remove the guidewire from the IntraClude device can result in diminished cardioplegia flow.

WARNING: Prior to initiating cardiopulmonary bypass, completely

expel air from the system. Failure to eliminate air from the system could result in an air embolus.

WARNING: Ensure the RHV is completely closed to prevent air entering the IntraClude device.

11. Attach the fully primed cardioplegia line and aortic root venting line to the appropriate arms of the Y-connector. Remove the detachable hose barbs as required (see Figure 2).

12. With the white tubing clamp closed, open the cardioplegia and vent lines and circulate cardioplegia solution until the Y-connector is de-aired and the cardioplegia line is fully primed.

WARNING: Always use a vacuum-relief valve in the aortic root vent line. Overly vigorous venting or not using a vacuum-relief valve may cause introduction of air into the IntraClude device. In addition, overly vigorous venting may cause migration of the balloon toward the aortic valve.

WARNING: Prior to initiating cardiopulmonary bypass, completely expel air from the system. Failure to eliminate air from the system could result in an air embolus.

13. Close the cardioplegia line tubing clamp, leaving the aortic root vent tubing clamp open.

14. Pass off the red pressure tubing to anesthesia and attach to a transducer capable of measuring pressures up to 200 mmHg. Flush as needed and zero the transducer. Open the red stopcock on the IntraClude device to read the aortic root pressure.

15. Pass off the blue pressure tubing to the perfusionist and attach to a transducer capable of measuring pressure of at least 500 mmHg. Flush as needed and zero the transducer. Leave the blue stopcock off the IntraClude device.

16. To prepare for bypass, ensure that all connections (stopcocks, pressure lines, etc.) are tight and not leaking, ensure that all lines have been properly de-aired, and zero the pressure monitoring systems.

17. Gradually initiate cardiopulmonary bypass while closely monitoring arterial line pressure. Continue to monitor the arterial line pressure throughout the case.

WARNING: When the balloon is advanced through the cannula the arterial inflow lumen is restricted. Immediately discontinue bypass at high arterial pressure to prevent vascular injury and/or tubing disruptions.

WARNING: High arterial line pressure may indicate an arterial dissection and/or placement of the arterial cannula in a false lumen. Failure to immediately discontinue bypass may propagate dissection and injure the vasculature and/or aorta.

18. Turn on the aortic root vent.

INFATING THE BALLOON

1. In preparation for balloon inflation, connect the syringe (filled with sterile physiologic solution or diluted contrast media) to the side port on the blue stopcock. Close the stopcock to the pressure monitoring system.

CAUTION: If using contrast medium ensure that it has been diluted with sterile physiologic solution at a 6:1 saline:contrast ratio or difficulty inflating/deflating the balloon may occur.

WARNING: The patient should be hemodynamically stable with the heart decompressed and well drained prior to inflation of the balloon. Failure to fully decompress the heart prior to balloon inflation may result in movement of the balloon during inflation.

2. Use TEE and/or fluoroscopy to ensure that the balloon is correctly placed in the aorta above the sinotubular junction.

WARNING: Ensure the balloon is correctly placed in the aorta to prevent damage to the aortic valve or inadequate cerebral perfusion.

3. Gently pull back on the IntraClude device to remove slack and to anchor it onto the lesser curve of the aorta. Under TEE and/or Fluoroscopy visualization, watch for balloon movement away from the aortic valve. The slack has been removed once a direct relationship between movement at the proximal hub and the balloon is seen.

WARNING: Avoid excessive shaft retraction, as vascular injury and/or inadequate cerebral perfusion could result if the brachiocephalic junction is blocked.

WARNING: Avoid excessive shaft retraction/traction as damage to the catheter may occur.

4. Under TEE and/or fluoroscopy guidance, inflate the balloon to create an occlusive seal. The balloon should be inflated per the following steps:

- a) Using the diameter of the aorta measured prior to the start of the procedure and the included compliance chart (see Figure 3), determine the estimated target balloon diameter and approximate inflation volume.

- b) Begin to inflate with approximately 10 mL of sterile physiologic solution while watching for occlusion. Note - occlusion is obtained when the aortic root pressure is zero and can be verified by observing for leakage around the balloon using TEE.

- c) Slowly inject an additional 2 mL of sterile physiologic solution; watch for occlusion.

- d) Continue the previous step until occlusion has occurred.

- e) Once occlusion has been obtained, close the balloon inflation (blue) stopcock to the syringe and open it to the balloon pressure monitoring system.

Note: A typical initial balloon pressure is between 300 to 400 mmHg. Balloon pressure should not be used as an indicator of occlusion, but only as a reference for balloon functionality.

WARNING: Hold the device shaft taut during balloon inflation to reduce the risk of proximal balloon migration and potential damage to the aortic valve.

WARNING: Do not exceed maximum inflation volume of 40 mL as balloon burst may occur.

WARNING: Over inflation of the balloon, based on the aorta diameter and compliance chart (see Figure 3), may cause trauma to the aorta. Under inflation may cause inadequate occlusion.

WARNING: Always monitor pressure via the right radial or right

brachial arterial line to ensure adequate perfusion of the head and neck vessels. Cerebral oximetry can also be used to monitor the perfusion flow to the head vessels.

WARNING: Do not discontinue bypass at this time as migration of the balloon may occur once bypass is reinitiated. If a reduction in pump flow is desired do not reduce arterial return more than 50% of full flow.

5. Immediately following inflation, discontinue venting the aortic root. If cardiopulmonary bypass was reduced, resume standard cardiopulmonary bypass flow once occlusion has been obtained. Gradually begin infusion of cardioplegia. A rise in root pressure is expected during cardioplegia delivery.

WARNING: Ensure that the cannula or introducer sheath rotating hemostasis valve is not overtightened and that there is not excessive retraction/traction on the device or cardioplegia delivery can be restricted.

6. If balloon repositioning is required partially deflate, reposition, and re-inflate the balloon using TEE and/or Fluoroscopy guidance.

7. If the catheter is resting against the inner wall of the Y-arm, open the cannula hemostasis valve slightly and advance the catheter shaft (no more than 0.5 cm) in order to reduce tension in the shaft.

8. Close the rotating hemostasis valve again, taking care not to overtighten. Secure the device in place following the steps below.

- a) Secure the IntraClude device (Figure 2) in place by locking the Clamp-Lock device, while on the 10.5 Fr extended relief only, against the cannula hemostasis valve. The Clamp-Lock device must be twisted to the 'O' position to slide and to the 'I' position to lock.

WARNING: Avoid excessive shaft retraction/traction, as vascular injury and/or inadequate cerebral perfusion could result if the brachiocephalic junction is blocked.

WARNING: Avoid excessive shaft retraction/traction as damage to the catheter or reduction in cardioplegia delivery may occur.

WARNING: Ensure that the balloon is above the sinotubular junction as inadequate coronary perfusion of cardioplegia could result if the sinotubular junction is blocked.

WARNING: To reduce the risk of device migration and damage to the aortic valve, the Clamp-Lock device must be locked and against the cannula/introducer hemostasis valve.

WARNING: To reduce the risk of device migration and damage to the aortic valve, use the locking device on the intended area only (extended strain relief for the Clamp-Lock).

9. Subsequent infusions of antegrade cardioplegia solution through the IntraClude device should be given to maintain heart arrest and provide myocardial protection.

WARNING: To ensure proper cardioplegia delivery, close aortic root vent connection prior to infusion of cardioplegia solution.

WARNING: Ensure that the aortic root pressure during cardioplegia delivery is less than the systemic arterial pressure or distal balloon migration may occur.

Note: Careful monitoring of the ECG will determine whether cardioplegia delivery is effective. Early return of cardiac activity suggests inadequate initial cardioplegia, inadequate aortic occlusion, or inadequate venous drainage. If initial delivery of cardioplegia is ineffective, balloon integrity and position should be confirmed with TEE and/or fluoroscopy.

10. Vent the aortic root, to remove fluid and air, as needed during the procedure.

WARNING: Overly vigorous venting may cause introduction of air into the IntraClude device and/or migration of the balloon towards the aortic valve.

WARNING: A loss of occlusion (evident by blood in the surgical field and/or a rise in aortic root pressure) along with a significant loss of pressure within the balloon may indicate compromised balloon integrity. Completely deflate the balloon and remove the IntraClude device carefully according to the withdrawal instructions to reduce the risk of patient injury.

WARNING: Do not add volume to the balloon based on a gradual balloon pressure drop as bursting of the balloon can occur. A gradual balloon pressure drop is normal, due to material relaxation of the balloon, and does not indicate loss of occlusion in the aorta.

DEFLATING AND WITHDRAWING THE INTRACLUDER DEVICE

Once the surgical procedure has been completed and the heart is ready to be resuscitated, complete the following steps to de-air the aortic root and to safely deflate the balloon and remove the IntraClude device.

1. Turn on the aortic root vent. Using TEE visualization, examine the aortic root for air.

2. Once all of the air has been removed, momentarily reduce bypass while keeping the aortic root vent on.

3. Connect an empty syringe to the blue stopcock and close the stopcock to the pressure line. Pull back on the syringe and fully deflate the balloon. Resume complete bypass.

WARNING: Ensure that balloon deflation is complete by verifying that the volume withdrawn is equivalent to the volume injected and that the balloon pressure is 0 mmHg or less. Note – the balloon inflation stopcock must be turned on to the pressure line to measure the balloon pressure. Failure to do so may result in difficulty withdrawing the device into the arterial cannula and potential vessel injury.

4. Continue to vent the aortic root to ensure all air has been removed from the aortic root.

WARNING: Always monitor pressure via the right radial or right brachial arterial line, as the deflated balloon may block the head vessels. Cerebral oximetry can also be used to monitor the perfusion flow to the head vessels.

5. Once cardiac function has been restored, withdraw the IntraClude device until

Released Date: 2018-04-23

Status = Released

the double line shaft markings indicate that the balloon is beyond the tip of the arterial cannula (See mark 2 on Figure 1).

6. Once bypass has been discontinued, continue withdrawing the IntraClude device until the deflated balloon rests within the Y arm of the cannula, against the hemostasis valve.

CAUTION: If difficulty is encountered retracting the IntraClude device balloon into the arterial cannula, ensure the balloon is fully deflated and carefully withdraw the IntraClude device into the cannula. Do not withdraw the IntraClude device while deflating the balloon. Failure to do so may result in difficulty withdrawing the device into the arterial cannula and potential vessel injury.

7. Clamp the non-reinforced section of the arterial cannula.

CAUTION: Before clamping the tubing attached to the arterial cannula hemostasis valve, ensure that the IntraClude device balloon is pulled back against the hemostasis valve. Failure to do so may result in inadvertent clamping of the IntraClude device tip within the arterial cannula.

8. Open the cannula hemostasis valve and remove the IntraClude device.

9. De-air the cannula by momentarily releasing the tubing clamp and allowing for slight back bleeding.

Tighten the cannula hemostasis valve seal to full closure.

Ensure that no air has entered the system.

WARNING: Prior to initiating cardiopulmonary bypass, completely expel air from system. Failure to eliminate air from the system could result in an air embolus.

10. Remove tubing clamp.

CAUTION: Dispose of used product in accordance with established hospital protocols for biohazards to minimize risk of exposure to blood borne pathogens.

WARNING: This device is designed, intended, and distributed for single use only. Do not re-sterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

STORAGE AND HANDLING

Store sterile packaged device in a cool, dry place until ready to use. Do not expose to organic solvents, ionizing radiation or ultraviolet light, or alcohol-based fluids.

DEVICE DISPOSAL

Used devices should be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

TECHNICAL ASSISTANCE

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837

Outside the U.S. and Canada (24 hours): 949.250.2222

In the UK: 0870 606 2040, Option 4

In Ireland: 01 8211012, Option 4

Edwards, Edwards Lifesciences, the stylized E logo, EndoReturn, IntraClude, and ThruPort are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

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Edwards

Manufacturer 
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

949-250-2500
800-424-3278

CE 0344

EC REP
Edwards Lifesciences Services GmbH
Edisonstr. 6
85716 Unterschleissheim
Germany

04/18
10006915003 A
© Copyright 2018, Edwards Lifesciences LLC
All rights reserved.

**Vendor: Do Not Print This Page.
For Internal Edwards Lifesciences Use Only.**

 Edwards	Title: IFU, ICF100, EN-TR	Page: 61 of 61
	Label Item/Revision: 10006915003/A	
	Visual Representation/Revision: DOC-0087526/A	
	Drawing/Specification: DOC-0018984	
	English Master/Revision: DOC-0053888/B	
		Refer to supporting material for signatures.

NOTE

1. ALL ART PRINTS 100% BLACK UNLESS OTHERWISE NOTED.

INKS

BLACK

GRAYSCALE 55% BLACK